




There is a new option available for low-grade rectocele and cystocele patients.

Releve® is a new vaginal splinting device intended as a first-line therapy for patients to use at home, when needed.

Releve® employs a soft, silicone balloon that is inserted into the vagina and inflated by the user prior to bowel movements and voiding.



Why patients love Releve®?

- Dignified solution and doesn't require continuous wear/use.
- Reduces straining during bowel movements and voiding.
- Helps patients feel completely empty after a bowel movement and voiding.
- Easy to use, discreet and portable.
- Releve is latex-free.




Releve® should be replaced yearly. Disposable covers are provided and required for use to protect against bacterial contamination.

A first line treatment that doesn't require you to wear it 24/7. Since most symptoms only appear during bowel movements, Releve is designed to be used only when you need it.




How it's used:

- The inflated balloon provides mild, pressurized support to the weakened vaginal walls to prevent the formation of rectocele and/or cystocele.
- A squeeze bulbs allows patients to inflate the balloon to their preferred size.
- The handle allows patients to hold the balloon in place during a bowel movement or voiding.
- When complete, the device is deflated by pressing the air release valve. The device is removed and cleaned for next use.
- Releve includes single-use nitrile covers to prevent bacterial contamination with multiple uses.



Uninflated Device
with Protective Cover



Uninflated Device
with Protective Cover

How patients can obtain Releve:

- Wholesale providers such as Pelvic Floor Clinics or Physical Therapists.
- Order online at www.relevesupport.com (MSRP \$95.00).
- Order by phone by calling (714) 832-6300

Developed by colorectal surgeons as a solution for low-grade prolapse patients that are non-surgical candidates.

Releve can be used as stand-alone therapy for obstructed defecation / voiding, or in combination with other therapies, including diet modification, biofeedback or pelvic floor physical therapy.

Surgin Inc.
Tel: (714) 832-6300
Email: info@surgin.com
www.relevesupport.com
relevesupport.com/patent

To learn more about Releve please visit relevesupport.com or scan the QR code:

